



Bone Marrow Donation

Welsh Bone Marrow Donor Registry (WBMDR)

Welsh Transplantation & Immunogenetics Laboratory, Welsh Blood Service, Pontyclun CF72 9WB, Wales, U.K.
Telephone: +44(0)1443 622177. Fax: +44(0)1443 622176. e-mail: wbmdr@wales.nhs.uk

WBM-107A Rev 05/04/2019

PRESCRIPTION FOR STIMULATED HUMAN PERIPHERAL BLOOD STEM CELL COLLECTION

(To be completed by the transplant center)

Patient name:	Patient ID number: (assigned by patient's registry)
Transplant center:	Patient ID number: (assigned by donor's registry)
Donor ID number:	GRID:
Donor registry: WBMDR	

PRE-COLLECTION PERIPHERAL BLOOD SAMPLES (maximum 50 mls):

	mls EDTA		mls ACD	Other, please specify:
	mls Heparin		mls no anticoagulant	
Samples to be shipped to: Name: Address:			Invoice(s) to be sent to: Name: Address:	
NOTE: This blood will be shipped at the time of the donor physical exam unless otherwise requested.			NOTE: All invoices associated with the blood sample procurement / shipment, donor work-up and stem cell collection should be sent to this address for payment (list only the requesting hub's address).	
Phone no:			Phone no:	
Fax no:			Fax no:	
Email:			Email:	

STIMULATED PBSC COLLECTION:

Required CD34 pos. cells / kg	X 10 ⁶ /kg
X recipient weight (kg)	kg
= total number of CD34 pos. cells	X 10 ⁶
+ CD34 pos. cells for quality testing	X 10 ⁶
= Total number of CD34 pos. cells	X 10⁶
IRB / Ethics Board (or equivalent) Approval	(Yes/No/Not Applicable or Date)

DISCLAIMER: The cell products collected from this donor are intended solely for the purpose of immediate therapeutic treatment for the above mentioned patient. Excess cells may be stored for future infusion for this patient. No other uses of these cells are permissible. Cells not used for the therapeutic treatment of the above mentioned patient must be disposed of properly. The donor center must be provided detailed information concerning the use and/or disposal of all portions of this cell product. By accepting these cells, the transplant physician also accepts these terms and conditions. Requests for deviations from these terms must be submitted in writing to the donor center for approval.

Preferred method of overnight storage (if needed) of apheresed product(s):	Donor Plasma Requirements: (Yes or No) If Yes, amount:
Transport Temperature: (Special packing materials such as gel packs must be provided by the transplant center unless alternative arrangements have been made with the donor or collection center) :	
Specification for transport of products for UK recipients: Products will be transported using a validated transport system. During transport the component surface temperature (CST) of each product shall be between 2-24°C for periods not exceeding 6 hours. If transported for a period exceeding 6 hours the CST will be maintained between 2-10°C within 12 hours and for the remainder of the journey up to 24 hours.	
Additional comments:	

PERIPHERAL BLOOD SAMPLES TO BE COLLECTED AT FIRST APHERESIS (maximum 50 mls)

	mls EDTA		mls ACD		mls Product Sample
	mls Heparin		mls no anticoagulant	Other:	
Additional comments: <@12@>					

Transplant physician:	Signature:	Date: (Day/Month/Year)
-----------------------	------------	---------------------------

