

WELSH TRANSPLANTATION & IMMUNOGENETICS LABORATORY**RENAL TRANSPLANT LIVING DONOR REQUEST FORM**Gwasanaeth Gwaed Cymru
Welsh Blood Service**PLEASE SEE INSTRUCTIONS ON REVERSE FOR COMPLETION OF FORM AND TRANSPORTATION OF SAMPLES TO THE LABORATORY**

Donor Details	
Surname:	Address: ADDRESSOGRAPH LABEL HERE Postcode:
Forename(s):	
Previous Name:	
Date of Birth :	
NHS No.:	
Hospital No.:	
Gender:	
Potential Recipient Details	
Surname:	Address: ADDRESSOGRAPH LABEL HERE Postcode:
Forename(s):	
Previous Name:	
Date of Birth :	
NHS No.:	
Relationship or donor to recipient:	
Tests Required (See reverse of form for blood sample requirements)	
Full Crossmatch <input type="checkbox"/>	ABO Titres <input type="checkbox"/>
Virtual Crossmatch <input type="checkbox"/>	ABO Group (Donor):
Match Grade <input type="checkbox"/>	ABO Group (Recipient): <i>If ABO groups of the patient and/or donor are unknown please provide additional samples for ABO titre testing and inform the laboratory if testing is required</i>
Tests Requested by :	Samples collected by:
Name (Please print):	Name (Please print):
Signature:	Signature:
Date:	Date : Time:
Report Required by (Please specify date):	
Laboratory Use only	TTID Donor: Recipient:
Date & Time Received:	

Download the current version of this form from: <https://portal.welsh-blood.org.uk/wtail/request-forms/>

WELSH TRANSPLANTATION AND IMMUNOGENETICS LABORATORY RENAL TRANSPLANT LIVING DONOR REQUEST FORM

SAMPLE REQUIREMENTS

Testing Required		Samples Required
Live Donor Testing	Full crossmatch	40mls EDTA 6mls EDTA-hand-labelled 6mls EDTA (<i>ABO incompatible only</i>)
	Virtual Crossmatch Match Grade only Altruistic Donor	10mls EDTA 6mls EDTA-hand-labelled
ABO Titres	Donor & Recipient	6mls EDTA

SAMPLE LABELLING / COMPLETION OF REQUEST FORM

1. A separate request form **MUST** be completed for each person.
2. Samples and forms must each be labelled with a minimum of 3 points of identification. The details provided on the sample and form must match.
3. Addressograph labels may be used on both samples and forms **UNLESS** specified otherwise (see sample requirements table above).
4. **SAMPLES WHICH ARE NOT LABELLED CORRECTLY WILL BE DISCARDED.**

SAMPLE STORAGE & TRANSPORTATION

1. Samples for HLA-Typing, cross-matching and urgent rejection monitoring **MUST** be pre-arranged with the laboratory.
2. Urgent samples **MUST** be clearly marked and transported to the laboratory with minimum delay.
3. Sample transport is the responsibility of the requestor. The requestor must assure themselves that the chosen mode of transport is effective and appropriate to maintain sample integrity.
4. Samples should be stored/transported at ambient temperature (approximately 22°C) and delivered to the laboratory in a timely manner, preferably within 24 hours of collection.
5. Samples must be packaged in accordance with current EU Carriage of Dangerous Goods regulations and IATA regulations.

Further information

Telephone	Email	Fax
01443 622186/179	WTAILRenal@wales.nhs.uk	01443 622310
WHTN 01797 2186/2179	AntibodyScreening@wales.nhs.uk	WHTN 01797 2310